

p38 MAPK (6F5) Mouse mAb


ABZ3878

Product Overview

Name	p38 MAPK (6F5) Mouse mAb
Accession(Primary)	Q16539
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (Q16539)

Gene Name	MAPK14
About Gene Symbol	MAPK14 (Mitogen-Activated Protein Kinase 14) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and JAK-STAT/cytokine signaling. It is an important research target in inflammatory disease and is widely studied to elucidate cellular mecha



Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma Tissue using p38 (6F5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.