

p38 MAPK (5A1) Mouse mAb


ABZ3879

Product Overview

Name	p38 MAPK (5A1) Mouse mAb
Accession(Primary)	Q16539
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q16539)

Gene Name	MAPK14
About Gene Symbol	MAPK14 (Mitogen-Activated Protein Kinase 14) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and JAK-STAT/cytokine signaling. It is an important research target in inflammatory disease and is widely studied to elucidate cellular mecha



Western blot analysis of p38 (5A1) in C2C12, PC-12, HepG2, COS7, C6, 3T3 and Hela lysates using p38 MAPK antibody.