

Smad2 (6H5) Mouse mAb


ABZ3881

Product Overview


Name	Smad2 (6H5) Mouse mAb
Accession(Primary)	Q15796
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 52 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (Q15796)


Gene Name	SMAD2
About Gene Symbol	SMAD2 (SMAD Family Member 2) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Smad2 in HeLa, A431, Jurkat and K562 lysates using Smad2 antibody.



Immunocytochemistry analysis of Smad2 (6H5) in HeLa using Smad2 antibody.



Western blot analysis of Smad2 (6H5) in 3T3
lysates using Smad2 antibody.