


IKB beta (5H3) Mouse mAb

ABZ3882


Product Overview

Name	IKB beta (5H3) Mouse mAb
Accession(Primary)	Q15653
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (Q15653)

Gene Name NFKBIB**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma using IKB beta (5H3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using IKB beta (5H3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.