

IKB beta (3H3) Mouse mAb

ABZ3883


Product Overview

Name	IKB beta (3H3) Mouse mAb
Accession(Primary)	Q15653
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100


Target information (Q15653)

Gene Name NFKB1B

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma using IKB beta (3H3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using IKB beta (3H3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.