

SAFB1 (2E8) Mouse mAb

ABZ3884


Product Overview

Name	SAFB1 (2E8) Mouse mAb
Accession(Primary)	Q15424
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 103 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; ChIP-1:20


Target information (Q15424)

Gene Name SAFB

About Gene Symbol



Western blot analysis of SAFB1 in A431, MDA-MB-468, Breast cancer, ZR751 and MCF-7 lysates using SAFB1 antibody.



Immunocytochemistry analysis of SAFB1 (2E8) in HeLa using SAFB1 antibody.