

RSK1 p90 (6B9) Mouse mAb


ABZ3885

Product Overview


Name	RSK1 p90 (6B9) Mouse mAb
Accession(Primary)	Q15418
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 83 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q15418)

Gene Name	RPS6KA1
About Gene Symbol	RPS6KA1 (Ribosomal Protein S6 Kinase A1) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied



Western blot analysis of p90RSK in COS7, 3T3, Jurkat, HeLa and A549 lysates using p90RSK antibody.



Immunoprecipitation analysis of RSK1 p90 (6B9) in HeLa lysates using p90RSK antibody.