

Beclin 1 (4F4) Mouse mAb


ABZ3890

Product Overview

Name	Beclin 1 (4F4) Mouse mAb
Accession(Primary)	Q14457
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat, Hamster
Calculated MW	Calculated MW: 52 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q14457)

Gene Name	BECN1
About Gene Symbol	BECN1 (Beclin 1) is an autophagy regulatory protein involved in autophagy and apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Beclin 1 in HeLa, 3T3, C6, COS7, CHO-K1, rat brain and Jurkat lysates using Beclin 1 antibody.