


PTK2B (4B4) Mouse mAb

ABZ3891


Product Overview

Name	PTK2B (4B4) Mouse mAb
Accession(Primary)	Q14289
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 116 kDa; Observed MW: 116 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (Q14289)

Gene Name PTK2B**About Gene Symbol**

Western blot analysis of Pyk2 in A549, Ramos, Jurkat, THP-1 and Raji lysates using Pyk2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Tonsil tissues using PYK2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.