

HDAC1 (3F5) Mouse mAb


ABZ3892

Product Overview


Name	HDAC1 (3F5) Mouse mAb
Accession(Primary)	Q13547
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (Q13547)

Gene Name	HDAC1
About Gene Symbol	HDAC1 (Histone Deacetylase) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma Tissue using HDAC1 (3F5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using HDAC1 (3F5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.