

Olig2 (7A1) Mouse mAb

ABZ3893


Product Overview

Name	Olig2 (7A1) Mouse mAb
Accession(Primary)	Q13516
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q13516)

Gene Name OLIG2

About Gene Symbol



Western blot analysis of OLIG2 in CHO-K1 lysates and CHO-K1 transfected by EGFPOLIG2 fragment lysates using OLIG2 antibody.

