

beta III Tubulin (3H2) Mouse mAb


ABZ3894

Product Overview


Name	beta III Tubulin (3H2) Mouse mAb
Accession(Primary)	Q13509
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q13509)

Gene Name	TUBB3
About Gene Symbol	TUBB3 (Tubulin Beta 3 Class III) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in fibrotic disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of beta III Tubulin (3H2) in various cell lines using β 3Tubulin (3H2) antibody.



Western blot analysis of beta III Tubulin (3H2) in mouse brain, rat brain lysates using beta III Tubulin (3H2) antibody.