


clAP2 (1F12) Mouse mAb

ABZ3895


Product Overview

Name	clAP2 (1F12) Mouse mAb
Accession(Primary)	Q13489
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q13489)

Gene Name BIRC3**About Gene Symbol**

Western blot analysis of clAP2 in Ramos, COS7, Raji and Daudi lysates using clAP2 antibody.



Western blot analysis of cIAP2 (1F12) in Ramos, COS7 and Raji lysates using cIAP2 antibody.