

XRCC4 (7C10) Mouse mAb


ABZ3896

Product Overview


Name	XRCC4 (7C10) Mouse mAb
Accession(Primary)	Q13426
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38-45 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (Q13426)


Gene Name	XRCC4
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded XRCC4 in Human lung cancer using XRCC4 (7C10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of XRCC4 (7C10) in XRCC4 in T47D, 3T3, COS7, C6 and HeLa lysates using XRCC4 (7C10) antibody



Immunoprecipitation analysis of XRCC4 in HeLa lysates using XRCC4 antibody.