

ATM (9F7) Mouse mAb


ABZ3898

Product Overview


Name	ATM (9F7) Mouse mAb
Accession(Primary)	Q13315
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (Q13315)

Gene Name	ATM
About Gene Symbol	ATM (ATM Serine/Threonine Kinase) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded Human Colon Carcinoma Tissue using ATM (9F7) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using ATM (9F7) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.