

HP1 gamma (5G10) Mouse mAb


ABZ3899

Product Overview


Name	HP1 gamma (5G10) Mouse mAb
Accession(Primary)	Q13185
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat, Hamster
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (Q13185)


Gene Name	CBX3
About Gene Symbol	




Western blot analysis of HP1gamma in HeLa, 3T3, C6, COS7 and CHO-K1 lysates using HP1 gamma antibody.




Immunocytochemistry analysis of HP1 gamma (5G10) in HeLa using HP1 gamma antibody.



Immunoprecipitation analysis of HP1 gamma (5G10) in HeLa lysates using HP1gamma antibody.



Immunohistochemistry analysis of paraffin-embedded Prostate Cancer using HP1 gamma antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of HP1 gamma in HeLa using HP1gamma (2001375G10) antibody (green) and Histone H3.1 (Phospho-Ser10) (310045) antibody (red). DAPI (blue).