

CHD3 (7F6) Mouse mAb


ABZ3901

Product Overview


Name	CHD3 (7F6) Mouse mAb
Accession(Primary)	Q12873
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 227 kDa; Observed MW: 260 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (Q12873)


Gene Name	CHD3
About Gene Symbol	CHD3 (Chromodomain Helicase DNA Binding Protein 3) is a transcription factor involved in adaptive immunity. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of CHD3 (Cterminus) in NTERA-2, 3T3, C6, K562 and HeLa lysates using CHD3 (Cterminus) antibody.



Immunocytochemistry analysis of CHD3 (7F6) in HeLa using CHD3 (Cterminus) antibody.



Immunoprecipitation analysis of CHD3 (7F6) in
Hela lysates using CHD3 antibody.