

BTK (5B12) Mouse mAb

ABZ3905


Product Overview

Name	BTK (5B12) Mouse mAb
Accession(Primary)	Q06187
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 76 kDa; Observed MW: 76 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20


Target information (Q06187)

Gene Name BTK

About Gene Symbol



Western blot analysis of BTK in Daudi, Jurkat(BTK negative), K562, Ramos and Raji lysates using BTK antibody.



Immunoprecipitation analysis of BTK (5B12) in Raji lysates using BTK antibody.