


GAD65 (6D5) Mouse mAb

ABZ3907

Product Overview

Name	GAD65 (6D5) Mouse mAb
Accession(Primary)	Q05329
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q05329)

Gene Name GAD2**About Gene Symbol**

Western blot analysis of GAD2 in rat
Hippocampus and rat Brain lysates using GAD2
antibody.

