

## 14-3-3 eta /YWHAH(5F2) Mouse mAb

ABZ3908

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | 14-3-3 eta /YWHAH(5F2) Mouse mAb           |
| <b>Accession(Primary)</b>   | Q04917                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 28 kDa; Observed MW: 28 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (Q04917)

**Gene Name** YWHAH

**About Gene Symbol**



