

Cytokeratin 17 (8C6) Mouse mAb


ABZ3909

Product Overview


Name	Cytokeratin 17 (8C6) Mouse mAb
Accession(Primary)	Q04695
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (Q04695)


Gene Name	KRT17
About Gene Symbol	KRT17 (Keratin 17) is a keratin cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Immunohistochemistry analysis of paraffin-embedded Human liver tissue using Cytokeratin 17 (8C6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunofluorescence analysis of Cytokeratin 17 (8C6) in Human breast tissue using CK17 antibody (red), and DAPI (blue).



Western blot analysis of Cytokeratin 17 (8C6) in HeLa, MCF-7, 293T lysates using Cytokeratin 17 (8C6) antibody



Immunoprecipitation analysis of Cytokeratin 17 in
HeLa lysates using Cytokeratin 17 (8C6) antibody