


PPAR delta (4G5) Mouse mAb

ABZ3913


Product Overview

Name	PPAR delta (4G5) Mouse mAb
Accession(Primary)	Q03181
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (Q03181)

Gene Name PPARD**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human Brain Tissue using PPAR delta (4G5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using PPAR delta (4G5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.