

E2F1 (4G8) Mouse mAb


ABZ3915

Product Overview


Name	E2F1 (4G8) Mouse mAb
Accession(Primary)	Q01094
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (Q01094)


Gene Name	E2F1
About Gene Symbol	E2F1 (E2F Transcription Factor 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of E2F1 in C6, Raw264.7 and HeLa lysates using E2F1 antibody.



Immunofluorescence analysis of E2F1 (4G8) in HeLa using E2F1 antibody, and DAPI (blue).



Immunoprecipitation analysis of E2F1 (4G8) in HeLa lysates using E2F1 antibody.