

hnRNP U/p120 (1B11) Mouse mAb

ABZ3916


Product Overview

Name	hnRNP U/p120 (1B11) Mouse mAb
Accession(Primary)	Q00839
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 91 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q00839)

Gene Name HNRNPU

About Gene Symbol



Western blot analysis of hnRNP U in K562, Jurkat and HeLa lysates using hnRNP U antibody.

