


TFIIB (2F6) Mouse mAb

ABZ3918

Product Overview

Name	TFIIB (2F6) Mouse mAb
Accession(Primary)	Q00403
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q00403)

Gene Name GTF2B**About Gene Symbol**

Western blot analysis of TFIIB in Jurkat, MCF-7, COS7, 3T3 and HeLa lysates using TFIIB antibody.

