

Cytochrome C (7C10) Mouse mAb

ABZ3919


Product Overview


Name	Cytochrome C (7C10) Mouse mAb
Accession(Primary)	P99999
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Chicken
Calculated MW	Calculated MW: 12 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P99999)


Gene Name CYCS

About Gene Symbol






Western blot analysis of Cytochrome C (7C10) in HeLa, T47D, mouse brain, K562, rat brain lysates using Cytochrome C (7C10) antibody.



Immunocytochemistry analysis of Cytochrome C (7C1) (green) in HeLa using Cytochrome C (7C1) antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded human tonsil tissue using Cytochrome C (7C10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.