

Smad2/3 (6F7) Mouse mAb


ABZ3921

Product Overview


Name	Smad2/3 (6F7) Mouse mAb
Accession(Primary)	P84022
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52,60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100

Target information (P84022)


Gene Name	SMAD3/SMAD2
About Gene Symbol	SMAD3 (SMAD Family Member 3) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded Human Liver tissue using Smad2/3 (6F7) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Smad2/3 (6F7) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Smad2/3 (6F7) in C6, MCF-7, Jurkat and 3T3 lysates using Smad2/3 (6F7) antibody