

Histone H3 (Zebrafish Specific) (2C4) Mouse mAb

ABZ3924


Product Overview

Name	Histone H3 (Zebrafish Specific) (2C4) Mouse mAb
Accession(Primary)	P68431
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Zebrafish
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:500-1:1000


Target information (P68431)

Gene Name	H3C1
About Gene Symbol	H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el

Zebrafish skeletal muscle



Western blot analysis of Histone H3 (Zebrafish Specific) in Zebrafish skeletal muscle lysates using Histone H3 (Zebrafish Specific) (2C4) antibody



Western blot analysis of Histone H3 (Zebrafish Specific) (2C4) in C6, 3T3, HeLa lysates using Histone H3 (Zebrafish Specific) (2C4) antibody.