

Histone H3 (10A10) Mouse mAb


ABZ3925

Product Overview


Name	Histone H3 (10A10) Mouse mAb
Accession(Primary)	P68431
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P68431)


Gene Name	H3C1
About Gene Symbol	H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el



Western blot analysis of Histone H3 in HeLa, Raw 264.7, rat Testis lysates using Histone H3 (10A10) antibody



Western blot analysis of Histone H3 (10A10) in C6, 3T3, HeLa lysates using Histone H3 (10A10) antibody.



Immunoprecipitation analysis of Histone H3 (10A10) in Hela lysates using Histone H3 (1A1) antibody.