

ACTA1 (2B10) Mouse mAb


ABZ3926

Product Overview


Name	ACTA1 (2B10) Mouse mAb
Accession(Primary)	P68133
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (P68133)

Gene Name	ACTA1
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded rat liver tissue using alpha Skeletal Muscle Actin (2B10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Western blot analysis of alpha Skeletal Muscle Actin (2B10) in T47D, 3T3, COS7, C6 and HeLa lysates using alpha Skeletal Muscle Actin (2B10) antibody