


PP1C alpha (1C11) Mouse mAb

ABZ3930


Product Overview

Name	PP1C alpha (1C11) Mouse mAb
Accession(Primary)	P62136
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (P62136)

Gene Name PPP1CA**About Gene Symbol**

Western blot analysis of PPP1A in A431, Jurkat and HeLa whole lysates using PPP1A antibody.



Immunocytochemistry analysis of PP1C alpha (1C11) in Hela using PPP1A antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PPP1A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.