

## PP1C alpha (1H1) Mouse mAb

ABZ3931


### Product Overview


|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | PP1C alpha (1H1) Mouse mAb                 |
| <b>Accession(Primary)</b>   | P62136                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG1                                       |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 38 kDa; Observed MW: 38 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IF-1:50-1:200             |

### Target information (P62136)

**Gene Name** PPP1CA

**About Gene Symbol**





Immunocytochemistry analysis of PP1C alpha in HeLa cells using PP1C alpha antibody.