

SUMO2/3 (3F7) Mouse mAb


ABZ3932

Product Overview

Name	SUMO2/3 (3F7) Mouse mAb
Accession(Primary)	P61956
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 11 kDa; Observed MW: Refer to figures
Recommended Dilution	WB-1:500-1:1000

Target information (P61956)

Gene Name	SUMO2/SUMO3
About Gene Symbol	



Western blot analysis of SUMO2/3 in HeLa lysates using SUMO2/3 antibody.