

PPP4C (2F11) Mouse mAb


ABZ3934

Product Overview


Name	PPP4C (2F11) Mouse mAb
Accession(Primary)	P60510
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (P60510)


Gene Name	PPP4C
About Gene Symbol	PPP4C (Protein Phosphatase 4 Catalytic Subunit) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Protein Phosphatase 4C in HeLa and Jurkat lysates using Protein Phosphatase 4C antibody.



Immunocytochemistry analysis of PPP4C in HeLa cells using PPP4C antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Protein Phosphatase 4C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.