

Wnt3 (4F11) Mouse mAb

ABZ3936

Product Overview

Name	Wnt3 (4F11) Mouse mAb
Accession(Primary)	P56703
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P56703)

Gene Name WNT3

About Gene Symbol



