

HDAC4 (4A3) Mouse mAb


ABZ3937

Product Overview


Name	HDAC4 (4A3) Mouse mAb
Accession(Primary)	P56524
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 119 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P56524)


Gene Name	HDAC4
About Gene Symbol	HDAC4 (Histone Deacetylase 4) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HDAC4 in Jurkat, HeLa and K562 lysates using HDAC4 antibody.



Immunoprecipitation analysis of HDAC4 (4A3) in HeLa lysates using HDAC4 antibody.



Western blot analysis of HDAC4 (4A3) in C2C12, 3T3, K562, Jurkat and Hela lysates using HDAC4 antibody.