

PP-T Mouse mAb


ABZ3941

Product Overview

Name	PP-T Mouse mAb
Accession(Primary)	P53041
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P53041)

Gene Name	PPP5C
About Gene Symbol	PPP5C (Protein Phosphatase 5 Catalytic Subunit) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PP5 in HeLa, Jurkat, 293T, A431 and MCF-7 lysates using PP5 antibody.