


Hexokinase II (1E8) Mouse mAb

ABZ3942

Product Overview

Name	Hexokinase II (1E8) Mouse mAb
Accession(Primary)	P52789
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Mouse, Rat, Cat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 102 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P52789)

Gene Name HK2**About Gene Symbol**

Western blot analysis of Hexokinase II in 3T3, COS7 and C6 lysates using Hexokinase II antibody.

