


Hexokinase II (3B9) Mouse mAb

ABZ3943

Product Overview

Name	Hexokinase II (3B9) Mouse mAb
Accession(Primary)	P52789
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 102 kDa; Observed MW: 102 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P52789)

Gene Name HK2**About Gene Symbol**

Western blot analysis of Hexokinase II in COS7, K562 and HeLa lysates using Hexokinase II antibody.

