

Eg5 (4H3) Mouse mAb


ABZ3944

Product Overview


Name	Eg5 (4H3) Mouse mAb
Accession(Primary)	P52732
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 119 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (P52732)


Gene Name	KIF11
About Gene Symbol	KIF11 (Kinesin Family Member 11) is a kinesin motor protein involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Eg5 in MCF-7, 293T, Jurkat and HeLa lysates using Eg5 antibody.



Immunocytochemistry analysis of Eg5 (4H3) in HeLa using Eg5 antibody.



Immunoprecipitation analysis of Eg5 (4H3) in Jurkat lysates using Eg5 antibody.