


MAP2K6 (8A11) Mouse mAb

ABZ3945


Product Overview

Name	MAP2K6 (8A11) Mouse mAb
Accession(Primary)	P52564
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P52564)

Gene Name MAP2K6**About Gene Symbol**

Western blot analysis of MAP2K6 (8A11) in PC-12, C6, Ramos, Jurkat and HeLa lysates using MKK3/6 antibody.



Immunoprecipitation analysis of MAP2K6 (8A11)
in HeLa lysates using MKK3/6 antibody.