

Cyclin H (1A7) Mouse mAb

ABZ3946

Product Overview

Name	Cyclin H (1A7) Mouse mAb
Accession(Primary)	P51946
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P51946)

Gene Name CCNH

About Gene Symbol



