


Aldehyde dehydrogenase 10 Mouse mAb

ABZ3947

Product Overview

Name	Aldehyde dehydrogenase 10 Mouse mAb
Accession(Primary)	P51648
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P51648)

Gene Name ALDH3A2**About Gene Symbol**

Western blot analysis of Aldehyde dehydrogenase 10 in A549, U251 and HepG2 lysates using Aldehyde dehydrogenase 10 antibody.

