

MCM2 (2B5) Mouse mAb


ABZ3954

Product Overview


Name	MCM2 (2B5) Mouse mAb
Accession(Primary)	P49736
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 102 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20; FC-1:50-1:100

Target information (P49736)


Gene Name	MCM2
About Gene Symbol	




Western blot analysis of MCM2 in Jurkat, K562, 3T3 and HeLa lysates using MCM2 antibody.




Immunoprecipitation analysis of MCM2 (2B5) in HeLa lysates using MCM2 antibody.




Immunofluorescence analysis of MCM2 (2B5) in HeLa cells using MCM2 antibody(green),and DAPI (blue).



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using MCM2 (2B5) antibody.High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded ovarian cancer tissues using MCM2 (2B5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded colon cancer tissues using MCM2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.