

## beta Arrestin 1 (6A1) Mouse mAb


ABZ3955

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | beta Arrestin 1 (6A1) Mouse mAb            |
| <b>Accession(Primary)</b>   | P49407                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG2b                                      |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 47 kDa; Observed MW: 51 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (P49407)

|                          |                                                                                                                                                                                                                                                                                                          |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | ARRB1                                                                                                                                                                                                                                                                                                    |
| <b>About Gene Symbol</b> | ARRB1 (Arrestin Beta 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of betaArrestin 1 in K562 and HeLa lysates using betaArrestin 1 antibody.