

PINCH (3C12) Mouse mAb


ABZ3958

Product Overview


Name	PINCH (3C12) Mouse mAb
Accession(Primary)	P48059
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; FC-1:50-1:100

Target information (P48059)


Gene Name	LIMS1
About Gene Symbol	




Immunoprecipitation analysis of PINCH in HeLa lysates using PINCH antibody.



Western blot analysis of PINCH (3C12) in C6, Jurkat, K562, 3T3 and HeLa lysates using PINCH (3C12) antibody



Immunofluorescence analysis of PINCH (3C12) in HepG2 cells using PINCH (3C12) antibody (green) and DAPI (blue).



Flow cytometry analysis of Hela stained with PINCH antibody (blue) and negative control (red).