

Ki67 (10B6) Mouse mAb

ABZ3960


Product Overview

Name	Ki67 (10B6) Mouse mAb
Accession(Primary)	P46013
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Recommended Dilution	IHC-1:50-1:100

Target information (P46013)

Gene Name MKI67

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human lymph node (A), esophagus (B), lung cancer (C), rectum cancer (D), using Ki67 (10B6) antibody with DAB staining. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.