

JNK1 (1A4) Mouse mAb


ABZ3961

Product Overview


Name	JNK1 (1A4) Mouse mAb
Accession(Primary)	P45983
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 46,54 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (P45983)


Gene Name	MAPK8
About Gene Symbol	MAPK8 (Mitogen-Activated Protein Kinase 8) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of JNK1 in C6 and 3T3 lysates using JNK1 antibody.



Immunofluorescence analysis of JNK1 (1A4) in 3T3 using JNK1 antibody.



Western blot analysis of JNK1 (1A4) in CHO-K1 cell lysates(B) and CHO-K1 transfected by JNK1 fragment fusion protein(A) cell lysates using JNK1 antibody.