


ZAP70 (4B2) Mouse mAb

ABZ3965


Product Overview

Name	ZAP70 (4B2) Mouse mAb
Accession(Primary)	P43403
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (P43403)

Gene Name ZAP70**About Gene Symbol**

Western blot analysis of ZAP70 in CEM and Jurkat lysates, negative in the Ramos lysates using ZAP70 antibody.



Immunoprecipitation analysis of ZAP70 (4B2) in Jurkat lysates (ZAP70 positive expression cell line) and Ramos lysates (ZAP70 negative expression cell line) using ZAP70 antibody.