

STAT5a (5C10) Mouse mAb


ABZ3968

Product Overview

Name	STAT5a (5C10) Mouse mAb
Accession(Primary)	P42229
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 91 kDa; Observed MW: 91 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P42229)

Gene Name	STAT5A
About Gene Symbol	STAT5A (Signal Transducer And Activator Of Transcription 5A) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Stat5a in 3T3, A549, Molt4, K562 and HeLa lysates using Stat5a antibody.